

AVICHE

AH-C10

Medical Compressor Nebulizer



Medical grade



Patent technology



Small particles



Double adjustment



3 years warranty



Quality compressor



Low residual




Low noise




Professional Medical Grade

1 MMAD of Spray particles is around 3 μ m , 85% effective particles make medicine directly reach bronchus and alveolus pulmonis and improve medicine absorption.



New Upgrade Of Mask

2 Fit the face, soft material is not irritating in



Dual Adjustment Technique

3 Adjust the atomization speed according to the doctor's advice, scientific speed and reasonable adjustment of the particle size of the drug spray according to the treatment site.



Aviation Compress Technique

4 High efficient nebulization, Nebulization rate is >0.25ml/min



LoW Noise Design Technology

5 Low noise design make noise is below 60dB, reduce noise disturbnace.



Designed For Hospital



Designed For Family

Product Name:	Medical Compressor Nebulizer
Product Model:	AH-C10
Spay Particle Diameter:	MMAD about 3 μ m
Nebulization Rate:	>0.25ml/min
Compressor Air Flow:	>8.5L/min
Liquid Bottle Volume:	7ml
Size (L x W x H):	178 x 109 x 230 mm
Net weight:	2.2 kgs
Noise:	<60dB