

AVICHE

Medical **A35** Oxygen Concentrator

A35 is a new generation oxygen concentrator for homecare and medical usage, which use our patented technology and have the ability to adjust the flow rate at 0.1L/min and moreover, there is a large display which can real-timely display oxygen concentration and oxygen flow. All the designs are customized to meet precise operation requirement and enhance the customer's experience.



Take care of family health with oxygen




Relentless Pursuit of Higher Quality

Through breaking the technological innovation,experiencing a rigorous assessment of monitoring. AVICHE finally realized successful upgradation in product appearance,design and technology.

1

Stable: High precision oil-free compressor


Adopted germany technology, high efficiency,low noise,low temperature,durable.



2

Oxygen concentration: 3L can reach over 90%

A35 adopted advanced aviation technology and use double molecular sieve to make oxygen, high oxygen concentration, steady running and durable



3

Release 10 million pcs/cm³ negative ions

A35 equipped negative IONs generator, release 10 million pcs/cm³ negative ions, Firstly purify air and then make oxygen to make users enjoy more healthy life

4

Large LCD screen

Large touch screen make it easy to use 8 inches LCD touch screen display makes font clear and accuracy-control .

5

Infrared wireless remote control

Users can control the oxygen concentrator at different status, even lie on a bed or far away from the oxygen concentrator



Product model	A35
Oxygen flow	1.0-5.0L/min
Oxygen concentration(3L)	>92%
Input power	220V/110V 50Hz/60Hz
Nebulizer	✓
Wattage	240VA
Weight	11.3 kg
Maximum noise level	<55dB
Size	290 x 293 x 588 mm